D I G I T

# Digitensio PRO

BP 3BD1-4D

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# Wrist Watch Blood Pressure Monitor

## Instruction Manual

**This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.**

Read the instructions carefully before using this device.

Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.

Keep dry.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Type BF applied part.

## **Table of contents**

### **1. Introduction**

---

- 1.1. Features of BP 3BD1-4D
- 1.2. Important information about self-measurement

### **2. Important information on the subject of blood-pressure and its measurement**

---

- 2.1. How does high/low blood-pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done, if regular high/low values are obtained?
- 2.4. MAM technology (Microlife Average Mode)
- 2.5. IHD - Pulse Arrhythmia Detection

### **3. The various components of the blood-pressure monitor**

---

### **4. Putting the blood-pressure monitor into operation**

---

- 4.1. Inserting the batteries
- 4.2. Sleep mode of the device
- 4.3. Set time and date
- 4.4. Select the user
- 4.5. Select Measuring Mode: Standard or Averaging Mode

### **5. Carrying out a measurement**

---

- 5.1. Before the measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measuring procedure
- 5.5. Memory storage and recall of the measurements
- 5.6. Discounting a measurements
- 5.7. Delete memories

### **6. Error messages/malfunctions**

---

### **7. Care and maintenance, recalibration**

---

### **8. Guarantee**

---

### **9. Technical specifications and international standards**

---

**Intended use:**

This oscillometric blood pressure monitor is intended for measuring non-invasive blood pressure in people aged 12 years or older.

**1. Introduction**

---

**1.1 Features of your BP 3BD1-4D**

- Your blood pressure monitor BP 3BD1-4D is a fully automatic, digital blood pressure measuring device for use on the wrist. This monitor enables very fast and reliable measurement of the systolic and diastolic blood pressure as well your pulse by use of oscillometric method.
- This device is easy to use and is proven in clinical studies to provide excellent accuracy. The large display makes it easy to visualize operational status and cuff pressure during a measurement.
- Please read through this instruction manual carefully before use and then keep it in a safe place. For further questions on the subject of blood pressure and its measurement, please contact your doctor.

**ATTENTION!****1.2 Important information about self-measurement**

- The measurement results given by this device is not a diagnosis. It is not replacing the need for the consultation of a physician, especially if not matching the patient's symptoms. Do not rely on the measurement result only, always consider other potentially occurring symptoms and the patient's feedback. Calling a doctor or an ambulance is advised if needed.
- Self-measurement means control, not diagnosis or treatment. Your values must always be discussed with your doctor. You should never alter the dosages of any medication without direction from your doctor.
- The pulse reading is not suitable for checking heart pacemakers!
- If you have been diagnosed with a severe Arrhythmia or irregular heartbeat, measurements made with this instrument should only be evaluated after consultation with your doctor.

**Electromagnetic interference**

Do not use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Keep a minimum distance of 3.3 m from such devices when using this device.

## **2. Important information on the subject of blood-pressure and its measurement**

### **2.1 How does high/low blood pressure arise?**

- Your blood pressure level is determined in the circulatory center of your brain. Your nervous system allows your body to adapt or alter blood pressure in response to different situations. Your body alters your pulse and the width of blood vessels through changes in muscles in the walls of blood vessels.
- Your blood pressure reading is highest when your heart pumps or ejects blood. This stage is called your Systolic Blood Pressure.
- Your blood pressure is lowest when the heart rests. (in-between beats) This is called your Diastolic Blood Pressure.

### **2.2 Which values are normal?**

- Blood pressure is too high if at rest, the diastolic pressure is above 90mmHg and/or the systolic blood pressure is over 140mmHg. If you obtain readings in this range, consult your doctor immediately. High blood pressure values over time damage blood vessels, vital organs such as the kidney and even your heart.
- Should the systolic blood-pressure values lie between 140mmHg and 160mmHg and/or the diastolic blood-pressure values lie between 90mmHg and 100mmHg, consult your doctor. Regular self-checks will be necessary.
- When blood-pressure values are too low, i.e. systolic values under 100mmHg and/or diastolic values under 60mmHg, consult your doctor.
- Even with normal blood-pressure values, a regular self-check with your blood-pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately.
- If you are undergoing medical treatment to control your blood pressure, please keep a record of the level of your blood pressure by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. **Never use the results of your measurements to independently alter drug doses prescribed by your doctor.**

Table for classifying blood-pressure values: World Health Organization (unit: mmHg)

| Colour      | Range                           | Systolic  | Diastolic | Recommendation                |
|-------------|---------------------------------|-----------|-----------|-------------------------------|
| Light green | blood pressure too low          | ↓ 100     | ↓ 60      | Consult your doctor           |
| Light green | blood pressure optimum          | 100 - 120 | 60 - 80   | self-check                    |
| Green       | blood pressure normal           | 120 - 130 | 80 - 85   | self-check                    |
| Yellow      | blood pressure slightly high    | 130 - 140 | 85 - 90   | Consult your doctor           |
| Orange      | blood pressure too high         | 140 - 160 | 90 - 100  | Seek medical advice           |
| Red         | blood pressure far too high     | 160 - 180 | 100 - 110 | Seek medical advice           |
| Dark red    | blood pressure dangerously high | 180 ↑     | 110 ↑     | Urgently seek medical advice! |

The higher value is the one that determines the evaluation. Example: a readout value between 150/85 or 120/98 mmHg indicates «blood pressure too high».

#### Further information

- If your values are mostly «normal» under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called «labile hypertension». In any case, please discuss the values with your doctor.
- **Correctly measured diastolic blood-pressure values above 120mmHg require immediate medical treatment.**

### 2.3. What can be done, if regular increased/low values are obtained?

- Please consult your doctor.
- Increased blood-pressure values (various forms of hypertension) are associated long- and medium term with considerable risks to health. This concerns the arterial blood vessels of your body, which are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can be the result. Furthermore, with long-term continuously increased blood-pressure values, the heart will become structurally damaged.
- There are many different causes of the appearance of high blood-pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organic malfunctions.

Please consult your doctor for information about the possible origins of your own increased blood-pressure values.

- d) There are measures which you can take, not only for reducing a medically established high blood-pressure, but also for prevention. These measures are part of your general way of life:

**A) Eating habits**

- Strive for a normal weight corresponding to your age. Reduce overweight!
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

**B) Previous illnesses**

Follow consistently any medical instructions for treating previous illness such as:

- Diabetes (Diabetes mellitus)
- Fat metabolism disorder
- Gout

**C) Habits**


- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (Coffee)

**D) Physical constitution**

- After a preliminary medical examination, do regular sport.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your sporting activities. He will advise you regarding the type and extent of types of sport that are possible for you.

## 2.4. MAM technology (Microlife Average Mode)

- «MAM technology» is a new type of concept for optimum reliability in self-measurement of blood-pressure.
- An advanced measurement accuracy is achieved by the automatic analysis of three successive measurements.
- The new system provides reliable values for the doctor and can be used as the basis for reliable diagnostics and medication therapy for high blood-pressure.

### A) Why MAM?

- Human blood-pressure is not stable

### B) Key advantages

Reduction of:

- Device scattering
- Insufficient rest prior to measurement
- Movement artefacts
- Cuff positioning influences

### C) Medical benefits

- Improved accuracy
- Reliable patient self-measurement data for the doctor
- Safe hypertension diagnostic
- Reliable therapy control

### D) Measurement sequence

- With the full measurement cycles, the total measurement time remains less than 3 minutes, compared to a single measurement time of 1.5 minute.
- Single results are not displayed.
- Due to the «Data Analysis» result, a 4th measurement may be applied.

## 2.5. IHD - Pulse Arrhythmia Detection

### Appearance of the Arrhythmia indicator

The appearance of this symbol  signifies that a certain pulse irregularity was detected during the measurement. The result can vary from your normal blood pressure. As a rule this is not a cause for concern; however, if the symbol appears more frequently (e.g. several times per week on measurements performed daily) or if it suddenly appears more often than usual, we recommend you inform your doctor.

Please show your doctor the following explanation:

#### Information for the doctor on frequent appearance of the Arrhythmia indicator


This BP 3BD1-4D wrist device is an oscillometric blood pressure measuring device that also analyses the pulse frequency during measurement. The accuracy of this device has been clinically validated.

If pulse irregularities occur during measurement, the arrhythmia symbol is displayed after the measurement. If the symbol appears **more frequently** (e.g. several times per week on measurements performed daily) **or if it suddenly appears more often than usual**, we recommend the patient to seek medical advice.

The device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

### 3. The various components of the blood-pressure monitor

---


The illustration shows the blood-pressure monitor BP 3BD1-4D, consisting of:  
Wrist cuff fits wrist size 13.5–19.5 cm

### 4. Using your blood pressure monitor

---

#### 4.1. Inserting the batteries

Insert batteries immediately after unpacking the unit. Battery compartment is located on the left lower side of the device (see illustration).

- Remove cover as illustrated:
- Insert the batteries (2 x size 1.5V), observing the indicated polarity
- If the battery warning  appears in the display, the batteries are empty and must be replaced by new ones. After battery warning  appears, the device will not work until batteries have been replaced.

## Attention!


- Please use «AAA» Long-Life or Alkaline 1.5V Batteries. Use of rechargeable batteries (1.2V Accumulators) is not recommended.
- If your blood-pressure monitor is left unused for long periods, please remove the batteries from the device.

**Functional check:** Press and hold the «M» button to test all the display elements. When functioning correctly all segments appear.

## 4.2. Sleep mode of the device

The device is designed to stay at the sleep mode after batteries are inserted. In this mode, the time and user icon will be displayed. You can conduct the following operations in this mode.


- Reading the set date: Press the «TIME» button. The date set will be displayed for 2 seconds.
- Recalling measurements: Press the «MEMORY» button. The measuring data stored can be displayed. For more information refer to the section on Memory storage and recall.


## 4.3. Setting time and date

This blood-pressure monitor incorporates an integrated clock with date display. This has the advantage, that at each measurement procedure, not only the blood-pressure values are stored, but also the exact moment of the measurement. After new batteries have been inserted, the device will request a user to set up time/date. (Example: Entering 2008-08-15 Time 09:30).

1. After new batteries have been inserted, the display comes into direct year setting, during which the 4 characters (2005) blink.

2. Enter the correct year by pressing «**MEMROY**» button. Each click will change the clock one year. Example: 3 x press.
3. Press «**TIME**» button. The display now switches to the current date, during which the first character (month) blinks.
4. Enter the current month by pressing the «**MEMORY**» button. Example: 7 x presses
5. Please follow the instructions above to set the day, hour and minutes.
6. After all settings have been made, click the «**TIME**» button once. The date is briefly displayed and then the time. The input is now confirmed and the clock begins to run.

## Note

- If you want to change date and time, press and hold the time button down for approx. 3 seconds till Year flashes, then you can enter the new values as described above.
- 1 x press on «**MEMORY**» button or «**TIME**» button advances one operation.
- Holding the button down speeds up the procedure switches to the current date, during which the first character (month) blinks.

## 4.4. Select the User

This blood pressure monitor is designed to store 30 measurements for each of two users. Before taking a measurement, be certain that the correct user has been selected.

- a) With the unit off, press and hold the «**TIME**» button for 3 seconds until the user icon in the lower left corner of the LCD screen flashes.

b) Press the «**MEMORY**» button to toggle between users.

c) Press the  button to make your selection.

#### 4.5. Select Measuring Mode: Standard or Averaging Mode

This instrument enables you to select either standard (standard single measurement) or Averaging Mode (automatic triple measurement).

To select standard mode, push the switch on the side of the instrument towards «**1**», then only one measurement will be taken. If you select «**3**», the unit will inflate and deflate three times resulting in one final average measurement.

#### Averaging Mode Technology

- In Averaging Mode, three measurements are automatically taken in succession and the result then automatically analyzed and displayed. Because your blood pressure constantly fluctuates, a result determined in this way is more reliable than one produced by a single measurement.
- By pushing mode switch to «**3**», the selected Averaging Mode appears in the display as the  symbol.
- The bottom, right hand section of the display shows a 1, 2 or 3 to indicate which of the 3 measurements is currently being taken.
- There is a break of 15 seconds between the measurements (15 seconds are adequate according to «Blood Pressure Monitoring, 2001, 6:145-147» for oscillometric instruments). A count down indicates the remaining time and a beep will sound 5 seconds before the 2nd and 3rd readings will begin.
- Do not remove the cuff between measurements.

## 5. Carrying out a measurement

---

### 5.1. Before measurement:

- Avoid eating, smoking as well as all forms of exertion directly before measurement. All these factors influence measurement result. Try to find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before measurement.
- Remove any garment that fits closely to your wrist arm.
- Always measure on the same arm (normally left).
- Take measurements regularly at the same time of day, since blood-pressure changes during the course of the day.


### 5.2. Common sources of error:

**Note:** Comparable blood-pressure measurements always require the same conditions! These are normally always quiet conditions.

- All efforts by the patient to support their arm can increase blood-pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during measurement. Use a cushion for support if necessary.
- If the wrist artery lies considerably lower (higher) than the heart, an erroneously higher (lower) blood pressure will be measured! Each 15cm difference in height results in a measurement error of 10mmHg!
- With repeated measurements, blood can accumulate in the respective arm, which can lead to false results. Correctly executed blood pressure measurements should therefore first be repeated after a 1 minute pause or after the arm has been held up in order to allow the accumulated blood to flow away. (Does not apply to a MAM measurement).

### 5.3. Fitting the cuff

- a) Remove all objects and jewelry (e.g. wristwatch) from the wrist in use. Put the cuff over the wrist with the monitor face on the inside of the wrist.
- b) The distance between the cuff and the hand should be approx. 10 mm (1/2").

- c) Secure the cuff with the Velcro fastener, so that it lies comfortably and not too tight, on the wrist.
- d) Lay your arm on a table, with the palm upwards. Support the arm a little with a cushion, so that the cuff rests at about the same height as the heart. Remain still for 2 minutes, before beginning the measurement.




### Note:

If it is not possible to fit the cuff to your left wrist, it can also be placed on the right. However, all measurements should be made using the same wrist.

### 5.4. Measuring procedure

After cuff has been appropriately positioned, the measurement can begin:

- a) Press the  button. The pump begins to inflate the cuff. Rising pressure in the cuff is shown in the display.
- b) After suitable inflation pressure has been reached, pump stops and pressure gradually falls, with cuff pressures being displayed. In case that inflation pressure is not sufficient, the monitor automatically re-inflates to a higher level.
- c) When the instrument detects a pulse, the heart symbol in the display starts to flash and a beep is heard for every heartbeat.
- d) A longer beep is sounded for 3 seconds when the measurement has been completed. The systolic, and diastolic blood pressures and pulse rate now appear in the display. Example (Fig.): Systole 129, Diastole 78, Pulse 75

- e) The appearance of this symbol  signifies that an irregular heart-beat was detected. This indicator is only a caution. It is important that you be relaxed, remain still and do not talk during measurements.

**NOTE:** We recommend contacting your physician if you see this indicator frequently.

- f) The measurement readings remain on the display until you switch off the instrument. If no button is pressed for a period of 1 minute the instrument switches itself off in order to preserve the batteries.

- g) When the unit is set to **Average Mode setting**, generally 3 separate measurements will take place in succession and calculates your detected blood pressure value. There will be 15 seconds resting time in-between each measurement. A count down indicates the remaining time and a beep will sound 5 seconds before the 2nd and 3rd readings will begin.

In case that the single data of each cycle differ too much from each other, a fourth measurement is performed before the result will be displayed. In rare cases the blood pressure is such unstable that even after four measurements the data vary too much. In this case «**ERR 6**» is shown and no result can be given. If one of measurement causes an error message it is repeated.

## 5.5. Memory – storage and recall of measurements

The blood-pressure monitor automatically stores the last 30 measurement values for each of 2 users. By pressing «**MEMORY**» button, the screen first shows the most recent readings along with time and date. Each time you press the «**MEMORY**» button an earlier measurement will be displayed along with time and date. If an irregular heartbeat interval was detected in any measurements stored in memory, the icon  will be displayed.  indicates measurement is taken under average mode.


## Further information

Measurements should not occur soon after each other, since otherwise, the results will be falsified. Wait therefore for several minutes in a relaxed position, sitting or lying, before you repeat a measurement.



### 5.6. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the  button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.


### 5.7. Delete memory

To delete stored memory readings, hold on «MEMORY» button down till LCD shows a symbol of «CL» flashing. Then press memory button again, you will hear 3 short beeps indicating memories have been deleted.

## 6. Error messages/malfunctions

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (Example: Error No. 2).

| Error No. | Possible cause(s)                                                                                                                                                                                                                     |
|-----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ERR 1     | No pulse has been detected.                                                                                                                                                                                                           |
| ERR 2     | Unnatural pressure impulses influence the measurement result. Reason: The arm was moved during the measurement (Artefact).                                                                                                            |
| ERR 3     | The inflation of the cuff takes too long. The cuff is not correctly seated.                                                                                                                                                           |
| ERR 5     | The measured results have indicated an unacceptable difference between the systolic and diastolic pressure. Take another measurement, carefully following the instructions. Contact your doctor, if you still obtain unusual results. |
| ERR 6     | Single data differ too much during Average Mode even after 4 cycles. No average result can be displayed.                                                                                                                              |

|           |                                                                                                                                                            |
|-----------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>HI</b> | The pressure in the cuff is too high (over 300mmHg) OR the pulse is too high (over 200 beats per minute). Relax for 5 minutes and repeat the measurement.* |
| <b>LO</b> | The pulse is too low (less than 40 beats per minute). Repeat the measurement.*                                                                             |

\* Please consult your doctor, if this or any other problem occurs repeatedly.

### Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

| <b>Malfunction</b>                                                                                                        | <b>Remedy</b>                                                                                                                                                                                                                                                                               |
|---------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The display remains empty when the instrument is switched on although the batteries are in place.                         | <ol style="list-style-type: none"> <li>1. Check whether the batteries are installed with the correct polarity and correct if necessary.</li> <li>2. If the display is unusual, re-insert batteries or exchange them.</li> </ol>                                                             |
| The device frequently fails to measure the blood pressure values, or the values measured are too low (too high).          | <ol style="list-style-type: none"> <li>1. Check the positioning of the cuff.</li> <li>2. Measure the blood pressure again in complete peace and quiet under observance of the details made under point 5.</li> </ol>                                                                        |
| Every measurement results in different values although the device functions normally and the values displayed are normal. | <ul style="list-style-type: none"> <li>• Please read the following information and the points listed under 5.2. «Common sources of error». Repeat the measurement. Please note: Blood pressure fluctuates continually so successive measurements will show some variability.</li> </ul>     |
| The blood pressure values measured differ from those measured by the doctor.                                              | <ul style="list-style-type: none"> <li>• Record the daily development of the values and consult your doctor about them. Please note: Individuals visiting their doctor frequently experience anxiety which can result in a higher reading than at home under resting conditions.</li> </ul> |

### Further Information


The level of blood-pressure is subject to fluctuations even with healthy people. Important thereby is, **that comparable measurements always require the same conditions (rest conditions)!**

You must consult your specialist dealer or chemist if there are technical problems with the blood-pressure instrument. **Never attempt to repair the instrument yourself!** Any unauthorised opening of the instrument invalidates all guarantee claims!

## 7. Care and maintenance, recalibration

---

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use petrol, thinners or similar solvent. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. The cuff must not be washed!
- d) Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.
- e) Never open the device! Otherwise the manufacturer calibration becomes invalid!

### Periodical recalibration

Sensitive measuring devices must from time to time be checked for accuracy. We therefore recommend a periodical inspection of the static pressure display **every 2 years**.

Your specialist dealer would be pleased to provide more extensive information about this.

## 8. Guarantee

---

This blood-pressure monitor BP 3BD1-4D is guaranteed for **5 years** from date of purchase. This guarantee includes the instrument and the cuff but does not apply to damage caused by improper handling, damage from leaking batteries, accidents, not following the operating instructions or alterations made to the instrument by third parties. The guarantee is only valid upon presentation of the guarantee card filled out by the dealer.

## 9. Technical specifications

---

|                                     |                                                                       |
|-------------------------------------|-----------------------------------------------------------------------|
| <b>Weight:</b>                      | 148 g (with batteries & cuff)                                         |
| <b>Size:</b>                        | 68 (L) x 60 (W) x 70mm (H)                                            |
| <b>Storage conditions:</b>          | -20 to +55°C                                                          |
| <b>Humidity:</b>                    | 15 to 95% relative humidity maximum                                   |
| <b>Operating conditions:</b>        | 10 to 40°C                                                            |
| <b>Display:</b>                     | LCD-Display (Liquid Crystal Display)                                  |
| <b>Measuring method:</b>            | Oscillometric                                                         |
| <b>Pressure sensor:</b>             | Capacitive                                                            |
| <b>Measuring range:</b>             | SYS/DIA: 20 to 280 mmHg<br>Pulse: 40 to 200 beats per minute          |
| <b>Cuff pressure display range:</b> | 0–299 mmHg                                                            |
| <b>Memory:</b>                      | Automatically stores the last 30 measurements for each of the 2 users |
| <b>Measuring resolution:</b>        | 1 mmHg                                                                |
| <b>Accuracy:</b>                    | Pressure within $\pm 3$ mmHg<br>Pulse $\pm 5\%$ of the reading        |
| <b>Power source:</b>                | 2 dry alkaline cells (batteries) size AAA, 1.5 V                      |
| <b>Battery lifetime :</b>           | approx. 460 measurements (using new batteries)                        |
| <b>Accessories:</b>                 | Storage box; IB etc                                                   |
| <b>Device standard:</b>             | EN1060-1/-3/-4<br>IEC 60601-1; IEC 60601-1-2 (EMC); IEC 60601-1-11    |
| <b>IP Class:</b>                    | IP20                                                                  |
| <b>Expected service life:</b>       | Device: 5 years or 10'000 measurements<br>Accessories: 2 years        |

This device complies with the requirements of the Medical Device Directive 93/42/EEC.  
Technical alterations reserved.

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**Fabriqué par :**  Microlife Corporation
**Mandataire :** Microlife AG